

Therapeutic Plasma Exchange (TPE) Spectra Optia[®] Apheresis System

Version 12 Procedure Training
(including single-needle option)

Operator's Manual Information

Spectra Optia Apheresis System

Intended Use

The Spectra Optia Apheresis System, a blood component separator, may be used to perform the following therapeutic apheresis, cell collection, and cell processing procedures*:

- Therapeutic plasma exchange
- Therapeutic plasma exchange with a secondary plasma device
- Red blood cell exchange, depletion, and depletion/exchange
- Mononuclear cell collection from the peripheral blood
- Granulocyte collection from the peripheral blood
- White blood cell depletion
 - WBC reduction for patients with leukocytosis at risk for leukostasis (USA)
- Platelet depletion
- Processing of harvested bone marrow

*Procedure availability varies by country.

Operator's Manual Information Continued

Contraindications for Use

- Leukocytapheresis is contraindicated in AML FAB M3 (APL) because of the accompanying disseminated intravascular coagulation. (Vahdat L, et al., "Early mortality and the retinoic acid syndrome in acute promyelocytic leukemia: impact of leukocytosis, low-dose chemotherapy, PMN/RAR-alpha isoform and CD13 expression in patients treated with all-trans retinoic acid." Blood 1994; 84: 3843-3849. Daver, et al., "Clinical characteristics and outcomes in patients with acute promyelocytic leukaemia and hyperleucocytosis." British Journal of Haematology 2015, 168, 646-653.)
- Other contraindications for the use of the Spectra Optia system are limited to those associated with the infusion of solutions and replacement fluids as required by the apheresis procedure, and those associated with all types of automated apheresis systems.

Possible Adverse Events of Apheresis Procedures Include:

- Anxiety, headache, light-headedness, digital and/or facial paresthesia, fever, chills, hematoma, hyperventilation, nausea and vomiting, syncope (fainting), urticaria, hypotension, allergic reactions, infection, hemolysis, thrombosis in patient and device, hypocalcemia, hypokalemia, thrombocytopenia, hypoalbuminemia, anemia, coagulopathy, fatigue, hypomagnesemia, hypogammaglobulinemia, adverse tissue reaction, device failure/disposable failure, air embolism, blood loss/anemia, electrical shock hazard, fluid imbalance, inadequate separation of blood components.

Reactions to Blood Products Transfused During Procedures

- Reactions to transfused blood products can include fever, circulatory overload, shock, allergic reactions, alloimmunization, transfusion-related acute lung injury (TRALI), and graft-versus-host disease (GVHD), as well as transmission of infectious diseases and bacteria. (Sources: Circular of Information for the Use of Human Blood and Blood Components, AABB, et al, ed., April, 2006; Guide to the preparation, use and quality assurance of blood components, 10th Edition, Council of Europe Publishing; Toy P et al., "Transfusion-Related Acute Lung Injury: Incidence and Risk Factors." Blood, 2012; 119: 1757-1767.)

Restricted to Prescription Use Only:

- Operators must be familiar with the system's operating instructions.
- Procedures must be performed by qualified medical personnel.

Learning Objectives

After completing this training you will be able to do the following regarding a TPE procedure using the Spectra Optia system:

- Discuss the principles of the procedure.
- Enter and discuss the data needed to perform the procedure.
- Discuss how the data you entered affects the procedure and the run targets.
- View and change the data on the run values screen.
- Make changes to the data on the data, run, and end run menu screens.
- Optimize the run to achieve the desired procedure outcomes.
- Troubleshoot issues that may arise.
- Describe using the single-needle option with a TPE procedure.
- Understand the issues related to pediatrics/low total blood volume (TBV) patients.

Presentation Overview


- Introduction
- Preparing to Perform the Procedure
- Monitoring the Run
- Completing the Run
- Making Changes
- Optimization
- Troubleshooting
- Single-Needle Access
- Low-TBV Patients

Introduction

- Exchange Set
- Basic Principles of TPE
- Connector


Exchange Set

1. Replace line
2. Remove bag
3. Channel
4. Connector
5. Cassette
6. AC and saline tubing
 - AC Correct Connect luer*
 - Saline spike
 - Sterile barrier filters
7. AC line
8. AC check valve
9. Colored clamps




*Correct Connect availability is dependent upon regulatory approval

Basic Principles of TPE




Basic Principles of TPE – Channel

1. Anticoagulated whole blood enters the channel.
2. Red blood cells (RBC) flow to the reservoir.
3. Plasma is pumped to the reservoir or to the remove bag.



Connector

1. Plasma port
2. RBC port



Questions?

Preparing to Perform the Procedure

- Configuration: Medication Infusion Notification
- Configuration – TPE Procedures
- Configuration – Blood Warmer
- Channel Loading
- Single-Needle
- Patient Data
- Fluid Data
 - Replacement Fluid
 - Fluid Balance
- Run Values
- Spiking the Replacement Fluid
- Patient Connection

Configuration: Procedure

Config		Data	Run	End Run
System	Procedure	Report	Network	
Height Units	Weight Units	Medication Infusion Notification	Custom Prime Recommendation (% TBV)	
cm	kg	Yes	10	
Pressure Alarm Limit (mmHg)		AC Container		
Inlet	Return	Notification	Volume (mL)	
-250	400	Yes	750	
14:14 12-18-2018	Confirm			

Configuration – TPE Procedures

The screenshot displays the configuration interface for TPE procedures. The main menu includes 'Config', 'Data', 'Run', and 'End Run'. Under 'Config', there are sub-menus for 'System', 'Procedure', 'Report', 'Network', 'TPE', and 'Blood Warmer'. The 'TPE' sub-menu is selected, showing the following parameters:

AC Infusion Rate (mL/min/L TBV)	Inlet:AC Ratio (__:1)	Plasma Volumes Exchanged	Custom Replacement Fluid (%)
0.8	10.0	1.0	0

At the bottom of the screen, the time and date are 9:25 and 8-14-2018. There is a 'Confirm' button, a back arrow, a 'TPE' label, and a 'No' symbol.

Configuration – Blood Warmer

Config Data Run End Run

System Procedure Report Network TPE Blood Warmer



Blood Warmer

Return Line Tubing Set (mL)

Yes 40

Replace Line



No

9:42 8-14-2018 Confirm   TPE

Channel Loading

Use the standard filler.

1. Centrifuge collar is in the correct position.
2. Notch on the locking pin is visible.
3. Optical reference is visible.
4. Channel sits flush with the groove.




Single-Needle Procedure

Single-Needle Connector


Male luer connection to patient

Female luer connection to inlet line



Female luer connection to return line

Converting to Single-Needle Access



Touch **Single Needle** on the options screen and follow the on-screen instructions.



Patient Data

Config Data Run End Run

9:33
8-14-2018

Confirm   TPE

Fluid Data – Replacement Fluid







Run Values


Config	Data	Run	End Run	
AC Infusion Rate	Inlet:AC Ratio (_:1)	Plasma Removed (mL)	Run Time (min)	Plasma Volumes Exchanged
1.0	10.0	2781	47	1.0
Flow Rate (mL/min)	AC	Inlet	Plasma	Replace
12.0	12.0	120.0	71.4	56.5
Current (mL)				
Target (mL)	558	5580	3299	2583
9:47 8-14-2018	Confirm		↩	⊘ TPE

Config	Data	Run	End Run	
AC Infusion Rate	Inlet:AC Ratio (_:1)	Plasma Removed (mL)	Run Time (min)	Plasma Volumes Exchanged
1.0	10.0	2781	70 ↑	1.0
Flow Rate (mL/min)	AC	Inlet	Plasma	Replace
8.0 ↓	8.0 ↓	80.0 ↓	47.6 ↓	37.7 ↓
Current (mL)				
Target (mL)	558	5580	3299	2583
9:48 8-14-2018	Confirm		↩	⊘ TPE

Spiking the Replacement Fluid

Config	Data	Run	End Run
1. Clamp second replace line.	4. Unclamp second replace line.	6. Reclamp second replace line.	
2. Spike saline/albumin.	5. Prime both lines.	7. Place line into detector.	
3. Squeeze drip chamber.			
			
14:27 12-27-2018	Continue		 TPE

Patient Connection



Questions?

Monitoring the Run

- Main Run
 - AIM Graphic
- Platelet Flush
- View Port

Main Run

Exchanging plasma.

2 min 63 min 65 min

Packing Factor 15

Current	AC	Inlet	Plasma	Replace
Flow Rate (mL/min)	8.0	80.0	51.5	42.2
Volume (mL)	15	151	60	33


Inlet -51 mmHg Return 120 mmHg

AC Infusion Rate 0.8 Inlet:AC Ratio 10.0

AC to Patient 4 mL


12:12 8-14-2018 TPE

Main Run – AIM Graphic




Platelet Flush

1. AIM system establishes the interface.
2. AIM system continuously monitors the interface position.
3. Buffy coat accumulates.



Platelet Flush (Continued)

4. System changes the pump flow rates to lower the position of the interface.
5. Buffy coat is returned to the patient.
6. Interface position is returned to normal after the platelet flush is completed.



View Port



Questions?

Completing the Run

- Run Targets Attained
- Rinseback and Disconnect
- Procedure Summary


Run Targets Attained

	Config	Data	Run	End Run
Run targets attained.				
		Target	Current	
Volumes Exchanged		0.9	0.9	
Run Time (min)		65	65	
Plasma Removed (mL)		2402	2402	
Replacement Fluid Used (mL)		2200	2200	
15:08 29-07-2019	Rinseback			TPE

Rinseback and Disconnect

Config | Data | **Run** | End Run

Rinseback is in progress.

Rinseback Time Remaining 5:05	Return 0 mmHg 	Rinseback Volume 11 mL
----------------------------------	---	---------------------------


Return Flow (mL/min)	54.2	Fluid Balance (mL)	-97
----------------------	------	--------------------	-----

12:36
19-10-2012

TPE

Config | Data | **Run** | End Run

1. Close inlet saline line.
2. Clamp inlet and return lines.
3. Disconnect patient lines.
4. Seal AC, saline, and bag lines.



13:23
8-14-2018



Unload

TPE

Procedure Summary

To calculate the patient's fluid balance, use the values on the procedure summary screen:

- +371 mL (AC Used)
- 2392 mL (Remove Bag)
- +1872 mL (Replacement Used)
- +3 mL (Tubing Set)
- +155 mL (Rinseback)
- 9 mL (Total)
- +0 mL Bolus (if given)
- 9 mL Patient's fluid balance



Questions?


Making Changes

- Configuration Menu
- Data Menu
- Run Menu
- End Run Menu

Data Menu

- Change Procedure
- Patient Data
- Fluid Data
- Alarm History
- Report

Patient Data



Fluid Data



Alarm History



Report

Config | **Data** | Run | End Run

Patient Data | Fluid Data | Alarm History | **Report**

Start Time | Procedure | Sent

	Start Time	Procedure	Sent
	Current		
✓	19-06-2019 13:26	RBCX	
✓	21-05-2019 07:21	TPE	
✓	24-01-2019 14:42	TPE	
✓	30-11-2018 11:07	CMNC	
✓	06-11-2018 10:52	TPE	

12:33
11-07-2019

Send [Refresh] TPE

Config | **Data** | Run | End Run

Patient Data | Fluid Data | Alarm History | **Report**

Start Time | Procedure | Sent

	Start Time	Procedure	Sent
	Current		
✓	19-06-2019 13:26	RBCX	
✓	21-05-2019 07:21	TPE	
✓	24-01-2019 14:42	TPE	
✓	30-11-2018 11:07	CMNC	
✓	06-11-2018 10:52	TPE	

12:37
11-07-2019

Send [Refresh] [No Access] TPE

Run Menu

- Main Run – Caution Status
- Exchange Status
- Operation Status
- Bolus
- Strobe
- Run Values
- Options

Main Run – Caution Status

Exchanging plasma.

7 min 51 min 58 min Packing Factor 15

Current	AC	Inlet	Plasma	Replace
Flow Rate (mL/min)	7.9	78.8	51.7	38.1
Volume (mL)	54	544	227	152

Inlet -9 mmHg Return 64 mmHg

AC Infusion Rate 1.1 Inlet:AC Ratio 10.0

AC to Patient 18 mL


12:32 11-07-2019 **Caution Status** TPE

Caution Status


Target or actual AC infusion rate exceeds 1.2 mL/min/L TBV.

12:13 2-21-2023 TPE

Exchange Status



Operation Status



Bolus

Config Data **Run** End Run

Exchange Status Operation Status **Bolus** Strobe Run Values Options

Saline

Spike saline for bolus.

Volume 50 mL


Flow Rate 60.0 mL/min

Current Bolus 0 mL

12:33
8-14-2018

Start Bolus ↩ TPE

Strobe



Run Values

Config	Data	Run	End Run
Exchange Status	Operation Status	Bolus	Strobe
Run Values	Options		
Maximum AC Infusion Rate	Inlet:AC Ratio (∞:1)	Plasma Removed (mL)	Run Time (min)
0.8	10.0	3123	67
Plasma Volumes Exchanged			
1.0			
	AC	Inlet	Plasma
Flow Rate (mL/min)	8.0	80.0	54.8
Current (mL)	103	1035	699
Target (mL)	536	5358	3640
			Replace
			45.4
			557
			2987
11:43 19-10-2012	Confirm	↩	TPE

Config	Data	Run	End Run
Exchange Status	Operation Status	Bolus	Strobe
Run Values	Options		
Maximum AC Infusion Rate	Inlet:AC Ratio (∞:1)	Plasma Removed (mL)	Run Time (min)
1.0	10.0	3123	48 ↓
Plasma Volumes Exchanged			
1.0			
	AC	Inlet	Plasma
Flow Rate (mL/min)	12.0 ↑	120.0 ↑	82.2 ↑
Current (mL)	114	1141	774
Target (mL)	536	5358	3640
			Replace
			68.3 ↑
			620
			2987
11:43 19-10-2012	Confirm	↩	TPE

Options

The screenshot displays the 'Options' menu within the 'Run' section of a medical device interface. The menu is organized into several rows, each with a title and one or more buttons. The 'Run' tab is highlighted in yellow at the top. The 'Options' button is also highlighted in yellow. The settings are as follows:




Setting	Value
Rinseback	Yes
Medication Infusion Notification	Yes
Custom Prime	No
Blood Warmer	No
Saline Rinse	No
Return Line	Yes
Tubing Set (mL)	40
Single Needle	No
Replace Line	No

At the bottom of the screen, the time and date are shown as '12:13 8-14-2018'. There is a 'Confirm' button, a back arrow button, and the text 'TPE'.

End Run Menu

- Rinseback, Disconnect, Run Targets

Rinseback, Disconnect, Run Targets




Questions?

Optimization

- Plasma Volumes Exchanged
- Fluid Balance
- AC to Patient


Plasma Volumes Exchanged

The number of plasma volumes exchanged determines the percentage of disease mediator removed.




Fluid Balance

Patient TBV 3000 mL, target fluid balance 100%



Fluid Balance

Patient TBV 3000 mL



AC to Patient

1.0 plasma volumes exchanged, AC infusion rate 0.8 mL/min/L TBV

TBV 1000 mL, 28% Hct

Inlet:AC ratio	10:1	15:1	20:1
AC used	120 mL	77 mL	57 mL
AC to patient	39 mL	36 mL	34 mL

TBV 3000 mL, 28% Hct

Inlet:AC ratio	10:1	15:1	20:1
AC used	359 mL	232 mL	171 mL
AC to patient	115 mL	107 mL	102 mL

TBV 5000 mL, 28% Hct

Inlet:AC ratio	10:1	15:1	20:1
AC used	599 mL	387 mL	285 mL
AC to patient	191 mL	178 mL	171 mL

Note: Terumo BCT does not recommend inlet:AC ratios above 15.

Questions?

Troubleshooting

- Inlet and Return Access Alarms
- Turbulence
- Semi-Automatic Mode
- High Interface
- Hemolysis
- Clumping

Inlet and Return Access Alarms

Inlet pressure was too low.

Pressure in inlet line was too low.

Alarm number: 3015 (12.0)

Inlet access was not properly positioned.

Inlet line was obstructed.

Inlet access was too small.

Continue

Rinseback

Disconnect

20:34
12-07-2017

Opt_alarm_Exchange_AlarmCode(3015)

TPE

Return pressure was too high.

Pressure in return line was too high.

Alarm number: 3034 (12.0)

Return access was not properly positioned.

Return line was obstructed.

Return access was too small.

Continue

Rinseback


Disconnect

20:35
12-07-2017

Opt_alarm_Exchange_AlarmCode(3034)

TPE

Turbulence



1. AIM image
2. Algorithm control icon

- Look through the view port to verify if turbulence is present in the connector.

Turbulence may be caused by:


- High inlet pump flow rate
- Low packing factor
- Platelet swirling
- Hyperviscosity/mild lipemia



- Decrease the inlet pump flow rate to increase the packing factor
- Do nothing
- Enter Semi-Automatic mode

Semi-Automatic Mode


- Semi-Automatic mode icon appears on the screen.
- AIM system no longer controls the interface position.



High Interface

“AIM system detected RBC interface near top of channel”

- Look through the view port to verify the position of the interface and verify the entered Hct.



Interface is near the top of the connector and the entered Hct is incorrect.




Lower the interface by increasing the entered Hct by 3% up to 3 times to avoid platelet loss.

Interface is near the top of the connector and the entered Hct is correct.





Touch **Retry** to resume the procedure.

Hemolysis



- Certain patient conditions may cause hemolysis.
- If hemolysis related to the patient's condition is suspected, verify the presence of hemolysis before disabling the RBC detector.

Clumping



If clumping is suspected:

- Decrease the inlet:AC ratio to 8:1.
- Process 100 mL of inlet volume.
- Verify clumping has resolved.
- Consider increasing the inlet: AC ratio to 10:1

Questions?


Single-Needle Access

- Single-Needle Connector
- Convert Access to Single-Needle
- Optimization

Single-Needle Connector

Male luer connection
to patient


Female luer
connection to inlet line




Female luer connection to
return line

Convert to Single-Needle Access

- With a blood warmer on the return line



- Without a blood warmer on the return line



Convert to Single-Needle Access

Exchanging plasma.

Exchange Status | Operation Status | Bolus | Strobe | Run Values | Options

7 (Circled Needle Icon) | Clock | Hourglass | 22 min | 37 min | 59 min | Packing Factor: 6

Current	AC	Inlet	Plasma	Replace
Flow Rate (mL/min)	12.0	120.0	78.2	58.6
Volume (mL)	166	1659	941	704



Inlet: -59 mmHg | Return: -5 mmHg | AC Infusion Rate: 0.8 | Inlet:AC Ratio: 10.0

AC to Patient: 28 mL

11:27 | 8-31-2018 | TPE



Optimization

- Inlet pump flow rate is set by the system.
- Procedure continues at a new inlet pump flow rate set by the system.



Optimization

- Inlet pump flow rate is set by the operator.
- System will continue the run at the same inlet pump flow rate.



Questions?

Low-TBV Patients

- Minimum Data Entry Limits
- AC Management
- Fluid Balance
- Custom Prime – RBC
- Custom Prime – RBC (60%)

Minimum Data Entry Limits

- Patient data
 - Height: 12 in or 30 cm
 - Weight: 5 lb or 2 kg
 - TBV: 300 mL
(The system will not calculate the TBV for weight < 25 kg.)
- Inlet pump flow rate
 - The allowable operator entered minimum inlet pump flow rate is 5 mL/min.
 - System may set up an inlet pump flow rate of less than 5 mL/min and the operator can confirm it.

AC Management

- AC infusion rate
 - AC infusion rate may need to be increased to achieve an inlet pump flow rate ≥ 5 mL/min.
- Inlet: AC ratio
 - Inlet: AC ratio needs to be kept at a value that maintains proper anticoagulation.

Configured AC infusion rate 0.8 mL/min/L TBV and inlet:AC ratio 10:1				
Patient TBV (mL)	300	400	500	600
Initial inlet pump flow rate	2.7	3.5	4.3	5.2
Increased AC infusion rate	1.5	1.2	1.0	0.9
Inlet pump flow rate	5.0	5.3	5.4	5.8


Fluid Balance

- Target fluid balance
 - Patient tolerance of the procedure
- Blood warmer
 - Patient comfort
- Custom prime
 - Improved tolerance of the volume of the extracorporeal circuit

Custom Prime – RBC

Config	Data	Run	End Run
Enter data for custom prime.			
<div style="border: 1px solid gray; padding: 5px; margin: 10px auto; width: 80%;"> <p style="text-align: center; background-color: #008080; color: white; padding: 5px; margin-bottom: 5px;">RBC</p> <p style="text-align: center; background-color: #008080; color: white; padding: 5px; margin-bottom: 5px;">Plasma</p> <p style="text-align: center; background-color: #008080; color: white; padding: 5px;">Albumin</p> </div>			
RBC Unit Hct (%)	Maximum Inlet Flow Rate (mL/min)	Volume (mL)	
<input type="text" value="55"/>	<input type="text" value="100"/>	<input type="text" value="230"/>	
11:24 19-10-2012	<input type="button" value="Confirm"/>	<input type="button" value="↩"/>	<input type="button" value="⊘"/> TPE

Custom Prime – RBC



Custom Prime – RBC (60%)

Patient		200 mL RBC	No blood prime	240 mL RBC	No blood prime
		No blood warmer		40 mL blood warmer	
TBV	Hct (%)	Change in patient Hct (%)			
300 mL	25	+5	-13	+8	-14
	30	+4	-14	+7	-16
	35	+3	-15	+6	-17
	40	+2	-16	+4	-18
600 mL	25	+3	-7	+5	-8
	30	+2	-8	+4	-9
	35	+2	-8	+3	-10
	40	+1	-9	+2	-11
1000 mL	25	+2	-5	+3	-5
	30	+2	-5	+2	-6
	35	+1	-5	+2	-6
	40	+1	-6	+2	-7

Note: The table indicates the approximate change in the patient's Hct immediately after custom prime.

Questions?

Thank You

TerumoBCT.com

