


Place cursor on number or click for description

Trima Accel® Auto RBC, Plasma Set

Catalog No. 80520



- | | | |
|-------------------------------|--------------------------------|-------------------------------|
| 1 Donor line | 15 Centrifuge pressure sensor | 29 Vent bag |
| 2 Sample bag | 16 Inlet line to centrifuge | 30 Air removal bag |
| 3 AC/draw/return manifold | 17 RBC line from centrifuge | 31 RBC bag |
| 4 Anticoagulant (AC) line | 18 Plasma line from centrifuge | 32 Auto RBC filter |
| 5 Sterile barrier filter | 19 Plasma pump header | 33 Auto RBC line |
| 6 Cassette | 20 Cassette label | 34 Frangible connector |
| 7 AC pump header | 21 Platelet pump header | 35 Replacement fluid line |
| 8 Draw/return pressure sensor | 22 Centrifuge loop | 36 Crossover line |
| 9 AC line | 23 Channel | 37 Crossover line check valve |
| 10 Return line | 24 Inlet port | 38 Sample bag manifold |
| 11 Draw line | 25 Collection chamber | 39 Channel clamp |
| 12 Inlet pump header | 26 Plasma collect line | 40 AC check valve |
| 13 Return reservoir | 27 Plasma bag | |
| 14 Return pump header | 28 Vent bag line | |

Trima Accel[®] Auto RBC, Plasma Set

Part Descriptions

1. **Donor line** – provides access to the donor for draw and return.
2. **Sample bag** – used for the collection of blood samples from the donor and the diversion of the first aliquot of blood.
3. **AC/draw/return manifold** – consists of the access to the injection site and the connections for the AC line (4), the draw line (11), and the return line (10).
4. **Anticoagulant (AC) line (with orange spike)** – carries AC from the AC bag to the cassette (6).
5. **Sterile barrier filter** – a 0.2-micron filter that prevents bacteria from entering the system, thereby maintaining a functionally closed environment for the collection of blood components.
6. **Cassette** – guides the flow of blood and products through the tubing set.
7. **AC pump header** – the tubing segment that fits into the AC pump.
8. **Draw/return pressure diaphragm** – allows the draw/return pressure sensor to monitor pressure at the donor access site.
9. **AC line** – carries AC from the cassette (6) to the AC/draw/return manifold (3).
10. **Return line** – carries blood components back to the donor.
11. **Draw line** – carries anticoagulated whole blood into the tubing set.
12. **Inlet pump header** – the tubing segment that fits into the inlet pump.
13. **Return reservoir** – holds uncollected components for return to the donor. Contains a return filter (200 micron) to prevent the return of microaggregates to the donor.
14. **Return pump header** – the tubing segment that fits into the return pump.
15. **Centrifuge pressure sensor** – detects high pressure in the centrifuge. Also monitors for high pressure in the Auto RBC filter during collection.
16. **Inlet line to centrifuge** – carries blood to the centrifuge.
17. **RBC line from centrifuge** – carries red blood cells from the centrifuge for collection or return to the donor.
18. **Plasma line from centrifuge** – carries plasma from the centrifuge for collection or return to the donor.
19. **Plasma pump header** – the tubing segment that fits into the plasma pump.
20. **Cassette label** – used by the Trima Accel system RBC detector to identify a tubing set as capable of collecting Platelet products.
21. **Platelet pump header** – the tubing segment that fits into the platelet pump.

22. **Centrifuge loop** – consists of the following:
 - Four-lumen tubing – carries fluid into and out of the channel.
 - Sleeves – used to reinforce the tubing at flex points.
 - Collars – used to secure the two ends of the loop in the centrifuge.
 - Bearings – the contact points between the centrifuge arm and the loop.
23. **Channel** – contains blood components as they are separated.
24. **Inlet port** – routes incoming blood into the channel.
25. **Collection chamber** – routes separated blood components to the appropriate collect lines.
26. **Plasma collect line** – carries the collected plasma to the plasma bag (27).
27. **Plasma bag** – 1 L bag that holds collected Plasma product.
28. **Vent bag line** – carries displaced air to and from the vent bag (29).
29. **Vent bag** – holds displaced air from the system.
30. **Air removal bag** – used to remove air from the RBC bag (31).
31. **RBC bag** – a bag for storage of concurrently collected RBC product.
32. **Auto RBC filter** – leukoreduces red blood cells during collection.
33. **Auto RBC line** – carries RBC storage solution to the RBC product post-collection and after the donor is disconnected.
34. **Frangible connector** – occludes the line to prevent flow of air or fluid; when broken, allows the flow of RBC storage solution into the RBC product bag(s).
35. **Replacement fluid line (with green spike)** – used to deliver replacement solution to the return reservoir that is then delivered to the donor during the procedure.
36. **Crossover line** – used to bypass the channel when adding RBC storage solution to the RBC product post-collection and after the donor is disconnected.
37. **Crossover line check valve** – works with the crossover line clamp to prevent blood from entering the crossover line during collection.
38. **Sample bag manifold** – consists of the access to the injection site and the connections for the draw line (11) and the sample bag (2).
39. **Channel line clamp** – isolates the channel during storage solution delivery to the RBC product.
40. **AC check valve** – works with the crossover line check valve to allow the system to monitor the pressure in the Auto RBC filter during RBC storage solution delivery.