


Spectra Optia® IDL Set

Catalog No. 12320



- 1 Inlet line
- 2 Diversion bag
- 3 Inlet line manifold
- 4 Inlet saline line clamp (red)
- 5 AC check valve
- 6 Return line
- 7 Return line manifold
- 8 Return saline line clamp (blue)
- 9 Anticoagulant (AC) line (orange)

- 10 Saline line (green)
- 11 Inlet line trap
- 12 Inlet pressure sensor diaphragm
- 13 Centrifuge pressure sensor diaphragm
- 14 Reservoir
- 15 Return pressure sensor diaphragm
- 16 Collect pressure sensor diaphragm
- 17 Centrifuge loop
- 18 Channel

- 19 Connector
- 20 Vent bag
- 21 Plasma line/Replace line
- 22 Spikes for replacement fluid (white)
- 23 Plasma bag
- 24 Collection bag
- 25 Collect line
- 26 Sample bulb assembly
- 27 Accessory line