


Spectra Optia® Exchange Set

Catalog No. 12220



- 1 Inlet line
- 2 Inlet line manifold
- 3 Inlet saline line clamp (red)
- 4 AC check valve
- 5 Return line
- 6 Return line manifold
- 7 Return saline line clamp (blue)
- 8 Anticoagulant (AC) line (orange)

- 9 Saline line (green)
- 10 Inlet line trap
- 11 Inlet pressure sensor diaphragm
- 12 Centrifuge pressure sensor diaphragm
- 13 Reservoir
- 14 Return pressure sensor diaphragm
- 15 Plasma pressure sensor diaphragm
- 16 Centrifuge loop

- 17 Channel
- 18 Connector
- 19 Vent bag
- 20 Remove bag
- 21 Remove line
- 22 Extra remove line
- 23 Replace line